

# About Your New TLSO Anterior Overlap Orthosis






## The Orthotic & Prosthetic Centers

Your new TLSO Anterior Overlap spinal brace is designed to provide you with secure, comfortable upper body support. The information below, along with the specific guidance from your orthotist, and/or physician, will tell you what to expect when using your spinal orthosis.

Follow your health care providers' instructions carefully to use this device safely and effectively; this document is designed to give you general "getting started" instructions. Please take the time to read this information carefully. Always contact your orthotist if you have any questions.

### Putting on your Anterior Overlap Orthosis

1. Before donning the orthosis make sure to wear a tight fitting cotton T-Shirt or stockinette. The material should be smooth; wrinkles under the orthosis can cause irritation or skin breakdown.
2. First identify top and bottom of the brace. Bottom of orthosis should be slightly contoured where brace fits over hips. While standing, wrap orthosis around waist, and align waist indentation. Feed straps through loops and attach loosely.
3. Lie down on bed, flex hips slightly by bending at the knees, loosen straps, and align orthosis with waist indentation between hip and ribs. Once the orthosis is positioned properly, tighten straps starting with middle and then alternate between top and bottom straps. Tighten straps as needed.



Please contact us if you have any questions.

#### **Boston**

260 Tremont Street  
Biewend Building, Level B  
Boston, MA 02116  
617 695 0101  
617 695 0222 / Fax

#### **Braintree**

197 Quincy Avenue - Route 53  
Braintree, MA 02184  
781 794 9991  
781 794 1769 / Fax

#### **Concord**

801 Main Street, Suite 8  
Concord, MA 01742  
978 371 1684  
978 371 7504 / Fax

#### **Methuen**

380 Merrimack Street, Suite 2A  
Methuen, MA 01844  
781 794 9991  
781 794 1769 / Fax

#### **Newton**

130 Rumford Avenue, Suite 106  
Newton, MA 02466  
617 916 1650  
617 916 2693 / Fax

#### **North Smithfield**

63 Eddie Dowling Highway  
North Smithfield, RI 02896  
781 794 9991  
781 794 1769 / Fax

#### **Plymouth**

124 Long Pond Road, Suite 3  
Plymouth, MA 02360  
774 773 9739  
774 773 9679 / Fax

#### **West Yarmouth**

126B MidTech Drive  
West Yarmouth, MA 02673  
508 775 2570  
508 775 7609 / Fax



**800 634 0606**  
oandpcenters.com

## Positioning Your Spinal Brace

Use the waist indentations on the spinal orthosis to position your brace. The waist should rest between the top of your hipbone and your lowest rib (waist). Be sure the waist on each side of the brace is properly located and the orthosis is centered at your midline.

## Proper Strap Tension

Straps should be tightened so they are comfortably snug. The object is to eliminate excess motion between your body and the orthosis, but NOT to overtighten the straps. Look for pre-marked lines your orthotist may have made on the straps for quick reference.

## Jacket Migration

Most spinal braces have a tendency to migrate upward. A slight change in the location of the spinal brace is expected when you sit up after having been lying down. Be sure the straps on the orthosis are secured firmly enough to avoid migration.

## Comfort Tips

Always wear your spinal brace over a cotton T-shirt that provides coverage in all of the contact areas of the brace. Try to keep the fabric smooth against skin and avoid fabric wrinkles or bunching. If necessary, change your T-shirt halfway through the day to increase your comfort.

## Areas of Redness

Some minor surface redness is common and should dissipate quickly. Check for localized redness after each wearing period. If redness or irritation persists call your orthotist for adjustments.

## Adjustments

Adjustments for fit and function can be made by appointment.

## Follow-up Care and Emergencies

Follow-up visits with your orthotist are vital for your healing – don't skip any scheduled appointments. Call our office with any questions or problems. For medical emergencies, call 911 or your local hospital.

Please contact us if you have any questions.

### **Boston**

260 Tremont Street  
Biewend Building, Level B  
Boston, MA 02116  
617 695 0101  
617 695 0222 / Fax

### **Braintree**

197 Quincy Avenue - Route 53  
Braintree, MA 02184  
781 794 9991  
781 794 1769 / Fax

### **Concord**

801 Main Street, Suite 8  
Concord, MA 01742  
978 371 1684  
978 371 7504 / Fax

### **Methuen**

380 Merrimack Street, Suite 2A  
Methuen, MA 01844  
781 794 9991  
781 794 1769 / Fax

### **Newton**

130 Rumford Avenue, Suite 106  
Newton, MA 02466  
617 916 1650  
617 916 2693 / Fax

### **North Smithfield**

63 Eddie Dowling Highway  
North Smithfield, RI 02896  
781 794 9991  
781 794 1769 / Fax

### **Plymouth**

124 Long Pond Road, Suite 3  
Plymouth, MA 02360  
774 773 9739  
774 773 9679 / Fax

### **West Yarmouth**

126B MidTech Drive  
West Yarmouth, MA 02673  
508 775 2570  
508 775 7609 / Fax



**800 634 0606**  
oandpcenters.com